

Fax

Date: May 14, 2024

Number of pages including cover sheet: 2

To: (929) 207-0142
Fax: +19292070142

From: Guelph Orthotics
Guelph Orthotics
undefined,
Guelph, ON,
CA N1L 0B7
Phone: +12262134654
Fax: +12262431085

Notes:
TEST FAX TESTING

Sent Via Ooma Office

Patient Verification Form

Important

Because of the unique features of the SPL, this orthosis has special instructions regarding usage, maintenance, subsequent follow up visits, and warranties.

Fillauer recommends this or a similar form be signed by the patient and kept on file.

Usage and Limitations of the Orthosis

Your orthosis is a custom orthopedic device especially fabricated and aligned for your needs. The orthosis will enable you to stand, ambulate, or to overcome other challenges attributable to disease or limb instability. It is essential that you the patient follow all verbal and written instructions which should be given to you by the orthotist regarding use, operation, and maintenance.

The orthosis should not be used in water or other corrosive environments and contact with acids or solvents should be avoided. It should not be exposed to temperatures above 122°F (50°C) since this may cause a change in the operation of the orthosis. The adjustment cable from the knee joint to the remote control should have some slack to prevent the plastic cable from cracking and breaking.

Proximal Satellite (Remote Control)

The 3 different modes of operation are controlled with a Proximal Remote Switch that has a slide lever and button. The automatic lock position is achieved when the slide lever is even with the top of the switch. When pulled up, the lock is immediately free for sitting and returns to the automatic lock position when released. When the button is pushed, the slide lever snaps into the lower position which locks the knee at all times.

Care and Maintenance

Frequently inspect the orthosis for wear and function. Your orthotist should be seen every 6 months for routine follow-up appointments.

Although the orthosis is intended for active ambulation you should use your best judgment to assess stability on uneven or inclined surfaces.

Any changes in skin condition or function of the orthosis should be communicated to the orthotist immediately. An interface sock such as the SmartKnit interface sock (with Coolmax® /Lycra® and antimicrobial fabric that inhibits bacteria, fungus, and odor) is recommended for limb protection and comfort.

Warranty

The warranty will be void if the item has been exposed to a corrosive environment, or if the item has been used in extremely abusive activities that could result in injury.

Patient Acknowledgement

I understand the Swing Phase Lock Orthosis instructions concerning the usage, operation, maintenance, subsequent follow up visits, and warranty as given to me by my orthotist.

Date

Patient Signature

WF106/11-02